


CE

Explore Genetic Technology to
Improve the Quality of Life



Pluslife

Product Catalog for Pet Diagnosis



Next-Generation
Molecular POCT



Pluslife Pet Diagnosis

< Rapid >

Onsite result in **7-35 mins**

< Accurate >

Comparable to **PCR**

< Simple >

1-2 mins hands-on time

< Powerful >

Up to **4 pathogens** in 1 test



Mini Dock

Rapid Molecular POCT System

PM003



Features

- Palm-sized, weighing less than iPhone 14 (~240g)
- Free from cold-chain
- Results displayed via Indicator/APP/PC
- Support LIS system
- Compatible with Smart Charging Station

Application Scenarios



Hospital/Clinic



Pet Shop




Testing Center


Test Procedure



Collect Sample



Mix Sample



Prepare Reaction Card



Run Test

Dock Pro 8

Eight-Channel POC Molecular Analyzer

PM008A



Features

- 8 test channels work independently
- Free from cold-chain
- Benchtop size and stackable
- 7" touchscreen with intuitive control
- Support LIS system
- Barcode reader and result printing

Application Scenarios



Hospital/Clinic



Pet Shop




Testing Center


Test Procedure




Collect Sample




Mix Sample



Prepare Reaction Card



Scan Sample



Run Test

Test Menu

Feline Infectious Diseases

- FHV-1/C.felis/M.felis
- Feline calicivirus(FCV) *
- Feline immunodeficiency virus/Feline leukemia virus (FIV/FeLV) *

Canine Infectious Diseases


- CDV/CPIV/CAV-2/Bb
- Canine coronavirus(CCoV) *
- Babesia/Ehrlichia/Anaplasma/Leishmania *

Feline/Canine Infectious Diseases

- Giardia *
- Parvovirus(CPV/FPV) *

* In development





FHV-1/C.felis/M.felis

Nucleic Acid Test Card




Specimen type: eye swab, nasal swab or oropharyngeal swab

LoD: 1000 copies/mL

Multi-pathogen detection: up to 3 pathogens in 1 test

Hands-on time: within 1-2 mins

A comparison study was performed between Pluslife and PCR



Source: Pet clinic (China)



CDV/CPIV/CAV-2/Bb

Nucleic Acid Test Card




Specimen type: eye swab, nasal swab or oropharyngeal swab

LoD: 1000 copies/mL

Multi-pathogen detection: up to 4 pathogens in 1 test

Hands-on time: within 1-2 mins

A comparison study was performed between Pluslife and PCR



Source: Pet clinic (China)




FCV

Nucleic Acid Test Card

Specimen type: eye swab, nasal swab or oropharyngeal swab

LoD: 1000 copies/mL

Hands-on time: within 1-2 mins



CPV/FPV

Nucleic Acid Test Card

Specimen type: rectal swab

LoD: 1000 copies/mL

Hands-on time: within 1-2 mins





Guangzhou Pluslife Biotech Co.,Ltd.

 +86-020-31703986

 market@pluslife.com

 www.pluslife.com

 Floor 3, Block E, Runhui Science & Technology Park, No.18 Shenzhou RD, Huangpu District,
Guangzhou, Guangdong, China



Website



LinkedIn



Line

PLPC-INT-20230907